


SUMMARY

- Next generation sequencing (NGS) can provide comprehensive and massively parallel analysis of molecular alterations in oncogenes and tumor suppressor genes associated with cancer pathway signaling.
- Complementary and scalable enrichment procedures for direct amplicon sequencing of “druggable” cancer genes were developed for FFPE and Fine Needle Aspiration (FNA) specimens across two orthogonal NGS platforms.
- PCR enrichment procedures enabled uniform coverage across SuraSeq™ 200, SuraSeq™ 500, SuraSeq™ 7500 and Ion AmpliSeq™ cancer gene panels, with NGS read depths of >1000X and detection of variants representing as few as 1-3% of reads.
- The results support the utility of high sensitivity, high resolution mutation assessments across thousands of loci in heterogeneous FFPE and FNA tumor specimens.

MATERIALS AND METHODS

Three FFPE and FNA compatible PCR-based enrichment panels were developed. The first two were multiplexed PCR assays that targeted 8 amplicons in 5 cancer genes or 35 amplicons in 16 cancer genes, including the most common mutations in the MAPK/ERK and PI3K/AKT pathways. The third included nearly 1000 amplicons from 52 cancer genes. Primers were designed to avoid known SNPs, repetitive sequences, and pseudogenes whenever possible, and included adaptor sequences to enable direct sequencing on either the Ion Torrent PGM or the Illumina GAIIx. FNA DNA inputs from 10 ng or FFPE DNA inputs from 10 ng to 2 ug were evaluated for PCR enrichment, and samples were barcoded up to 36/lane (GAIIx) or 13/chip (PGM). Workflows for GAIIx NGS required ~3-6 weeks, whereas sample processing on the PGM required <3 days.



SuraSeq™ 200

BRAF	NRAS
HRAS	PIK3CA
KRAS	

Table 1. A 5 gene, 8 amplicon panel represents 211 unique mutations indexed in the COSMIC database.

PGM	Orthogonal Confirmation	
	POS	Wt
POS	9	0
Wt	0	71

Figure 1. Ion Torrent PGM of 10 FFPE specimens demonstrates 100% concordance with orthogonal assays.




Figure 2. Total mapped reads on a 316 chip across 13 barcoded specimens using the Ion Torrent PGM. NTC = Non template control.

SuraSeq™ 500

ABL1	FGFR1	HRAS	MET
AKT1	FGFR3	JAK2	NRAS
BRAF	FLT3	KIT	PDGFRA
EGFR	RET	KRAS	PIK3CA

Table 2. A 16 gene, 35 amplicon panel represents 540 unique mutations and >95% of all mutations in these genes listed in COSMIC.

PGM	Orthogonal Confirmation	
	POS	Wt
POS	13	1*
Wt	1**	80

*Below LOD by Sanger sequencing. **Only detected with a probe-based assay with < 1% LOD. No biopsy was negative.

Figure 3. Ion Torrent PGM NGS of 20 FNA specimens demonstrates 98% concordance with confirmation assays, including Sanger sequencing.

Category	Intact DNA (N=27)	FFPE DNA (N=39)
Average Depth	38,500	24,500
Maximum	72,738	44,747
Minimum	1,784	3,432
Range within 5-fold of average	95%	95%
Median Variant	0.27%	0.32%

Table 4. Both intact cell line DNA and FFPE DNA demonstrate low levels of “background” base substitution variants (0.3%).

Ion AmpliSeq™

ABL1	ATM	CSF1R	ERBB4	FGFR3	HRAS	KDR	MLH1	NRAS	PTPN 11	SMARCB1	TP53
AKT1	BRAF	CTNNB1	FBXW7	FLT3	IDH1	KIT	MPL	PDGFRA	RB1	SMO	VHL
AKL	CDH1	EGFR	FGFR1	GNAS	JAK2	KRAS	NOTCH1	PIK3CA	RET	SRC	
APC	CDKN2A	ERBB2	FGFR2	HNF1A	JAK3	MET	NPM1	PTEN	SMAD4	STK11	

Table 4. Ion Torrent AmpliSeq™ is a 46 gene, 190 amplicon panel that represents 739 mutations.




Figure 5. Mutation frequencies in pooled cell line DNA (Table 5) are highly correlated across operators and sequencing chips. AmpliSeq™ multi-sample barcoding was not available for the samples processed.

SuraSeq™ 7500

ABL1	DNMT3A	GNAQ	MET	PTCH1	TP53
AKT1	EGFR	HIF1A	MPL	PTEN	VHL
AKT2	ERBB2	HRAS	NF2	PTPN 11	
BRAF	FES	IDH1	NOTCH1	RB1	
CDH1	FGFR1	IDH2	NPM1	RET	
CDK4	FGFR3	IKKBK	NRAS	SMAD4	
CDKN2A	FLT3	JAK2	PAX5	SMARCB1	
CEBPA	FOXL2	KIT	PDGFRA	SMO	
CREBBP	GATA1	KRAS	PIK3CA	SRC	
CTNNB1	GNA11	MEN1	PIK3R1	STK11	

Figure 6. A 52 gene, 981 amplicon panel represents 7500 unique mutations. Amplicons designed for genes shown in bold were sequenced across all coding exons.




Figure 8. Twelve known cancer gene mutations from 12 residual clinical FFPE samples are highly correlated after enrichment with SuraSeq™ 7500 and SuraSeq™ 500 panels and NGS on different sequencing systems.

PGM	SuraSeq™ 200	GAIIx				
		SuraSeq™ 500	Ion AmpliSeq™	SuraSeq™ 500	SuraSeq™ 7500 (100ng)	SuraSeq™ 7500 (1000ng)
SuraSeq™ 200	1	0.96	0.94	0.98	0.81	0.97
SuraSeq™ 500	0.96	1	0.90	0.96	0.80	0.94
Ion AmpliSeq™	0.94	0.90	1	0.96	0.93	0.95
SuraSeq™ 500	0.98	0.96	0.96	1	0.83	0.97
SuraSeq™ 7500 (100ng)	0.81	0.80	0.93	0.83	1	0.90
SuraSeq™ 7500 (1000ng)	0.97	0.94	0.95	0.97	0.90	1
Known Input	0.93	0.97	0.93	0.91	0.88	0.89

Figure 10. Mutation quantification using distinct PCR enrichment panels and NGS platforms is highly correlated using known input cancer gene mutations for pooled cell line DNA (Tables 3, 5).

CONCLUSIONS

- Three distinct PCR workflows enabled high depth enrichment of cancer-associated gene regions in FFPE and FNA DNA from residual clinical specimens.
- Mutation loads as low as 1-3% can be accurately identified in both cancer cell line and FFPE tumor DNA; “background” variant detection was only ~0.3%.
- PCR-based enrichment of cancer gene “hotspots” in FFPE and FNA tumors revealed concordance in mutation detection using orthogonal methods.
- Ion Torrent NGS successfully confirmed novel mutations from screening studies using the Illumina GAIIx, suggesting utility for high sensitivity orthogonal mutation confirmation using a second NGS system.
- SuraSeq™ cancer gene panels supported a streamlined protocol, low DNA inputs, multiplex target amplification, and, importantly, efficient multi-sample barcoding, even on the Ion Torrent PGM. Thus, SuraSeq™ focused gene panels offered high depth sequencing for multiple samples per run, compared to lower depth reads for singleplex AmpliSeq™.
- The three proposed enrichment approaches can accommodate both large-scale, whole exon mutation assessments in ~96 samples per run, as well as “hotspot” mutation analyses across 15-50 genes with a rapid turnaround time (<1 week).




Figure 7. SuraSeq™ 7500 recovers the full range of expected SNP fractions from pooled HapMap DNA in a Latin square titration, identifying 209 SNPs (for each of six cell lines) across 1254 data points. Mixing ratios were designed such that the majority of SNPs were <20% abundance. Expected SNP fraction is based on results from individually sequencing each HapMap cell line.




Figure 9. As little as 250 ng FFPE DNA supports high depth NGS following SuraSeq™ 7500 enrichment.